



Nail Cleaning  
Laser Device



User Instruction

Thanks for using this nail cleaning laser device. Hope it will take you good experience.

Operation:

- 1.Get the device out of the box and place it on the nail and press the button to turn on the device to work.
- 2.The green light will begin to blink and turn off soon during the device working time. The green light will blink 3 times when the time is up, and the device will turn off automatically.
- 3.The red light will blink 5 times when the battery is low.
- 4. It should be used on a daily basis
- 5. It should be used at least one time and not over 3 times per day.

**Note:** The device only works for 7 minutes per time and turn off automatically. It may not be used while charging.

Warning:

- 1.The device should be used in strict compliance with safety and operating precautions.
- 2.No side effects.
- 3.Direct or reflected laser radiation should be prevented from direct eye contact.
- 4.DO NOT watch the laser directly when the device is working.
- 5.Pregnant women should not use the device closed to womb.
- 6.Keep the device away from children.
- 7.Children under 18 years old should be treated under the supervision of an adult.

Maintenance:

The exterior side of the device may be cleaned with moistened cloth. Once in a few times of use, it is recommended to clean the rubber pads with alcohol. It should be cleaned before other people to use it.

Storage temperature: -20℃~50℃

Troubles shooting

Light blinks	Low battery warning	Charge the battery with USB cable with a 5V adaptor
Device does not work	Technical fault	Contact dealer

Technical Data:



Laser wavelength:905nm  
Blue light wavelength:470nm  
Laser intensity:25W  
Maximum pulse intensity:22W  
Charging Voltage:5V  
Charging time:3-4 hours